



Signex

Cannulated Compression Screw



OSTEONIC
The Medical Explorer

INTRODUCTION

Overview

Signex Cannulated Compression Screw system offers the two kinds(2.2, 3.0) of screw head size to treat variety indication.

Continuous discussion with surgeons and personnel working in the OR and solution-finding efforts thereafter have enabled us to develop the Signex Cannulated Compression Screw.

With this innovative product, we have been able to fulfill all of our valued customers' demands, ensuring the maximum quality and functionality level of the products we have to offer.

INDICATION

Fixation of fracture, osteotomies, fusions and arthrodeses

All In One Tray

D9Q.01020
Tray Only



INSTRUMENT

I-Handle



C8P.01010

Screw Driver Shaft

Ø2.2 SCREW(75mm)



D8P.22010

Ø3.0 SCREW(75mm)



D8P.22020

Drill Guide

Ø2.2 / 3.0 SCREW



D8P.02010

K-Wire

Ø2.2 SCREW (Ø0.8/100mm)



D8P.04080

Ø3.0 SCREW (Ø1.1/100mm)



D8P.04110

K-Wire Depth Gauge

Ø2.2 SCREW (Ø0.8 K-Wire)



D8P.03090

Ø3.0 SCREW (Ø1.1 K-Wire)



D8P.03120

Drill Bit

Ø2.2 SCREW (Ø1.8 / 115mm)



D8P.31170

Ø3.0 SCREW (Ø2.1 / 95mm)



D8P.31210

Screws system

- ▶ Two head sizes accommodate patient anatomical condition.
- ▶ Wide range of screw lengths covers variety indication.
- ▶ High precise screw thread from accurate manufacturing process.
- ▶ Screws are made of most biocompatible material of titanium.
- ▶ Star shaped head design secures stable insertion with enough strength.



Length	Ø2.2				Ø3.0			
	Thread Length	Part No.	Thread Length	Part No.	Thread Length	Part No.	Thread Length	Part No.
	Headless Short		Headless Long		Headless Short		Headless Long	
10	4	D1D.22100			4	D2D.30100		
11	5	D1D.22110			4	D2D.30110		
12	5	D1D.22120			4	D2D.30120		
13	5	D1D.22130			4	D2D.30130		
14	5	D1D.22140			4	D2D.30140		
15	5	D1D.22150			4	D2D.30150		
16	5	D1D.22160			4	D2D.30160	5	D2E.30160
17	5	D1D.22170			4	D2D.30170	6	D2E.30170
18	5	D1D.22180			4	D2D.30180	6	D2E.30180
19	5	D1D.22190			4	D2D.30190	7	D2E.30190
20	5	D1D.22200			4	D2D.30200	7	D2E.30200
21	5	D1D.22210			4	D2D.30210	8	D2E.30210
22	5	D1D.22220	8	D1E.22220	4	D2D.30220	8	D2E.30220
23	5	D1D.22230			4	D2D.30230	8	D2E.30230
24	6	D1D.22240	8	D1E.22240	5	D2D.30240	8	D2E.30240
25	6	D1D.22250			5	D2D.30250	8	D2E.30250
26	6	D1D.22260	8	D1E.22260	5	D2D.30260	10	D2E.30260
27	6	D1D.22270			6	D2D.30270	10	D2E.30270
28	6	D1D.22280	9	D1E.22280	6	D2D.30280	10	D2E.30280
29	6	D1D.22290			6	D2D.30290	10	D2E.30290
30	6	D1D.22300	10	D1E.22300	7	D2D.30300	12	D2E.30300
32			11	D1E.22320	7	D2D.30320	12	D2E.30320
34			12	D1E.22340	8	D2D.30340	14	D2E.30340
36			13	D1E.22360	9	D2D.30360	14	D2E.30360
38			14	D1E.22380	9	D2D.30380	16	D2E.30380
40			15	D1E.22400	10	D2D.30400	16	D2E.30400
45					15	D2D.30450	21	D2E.30450
50					20	D2D.30500	26	D2E.30500



OSTEONIC Co., Ltd.

1206Ho, 38, Digital-ro 29-gil, Guro-gu, Seoul, Korea P.O. #08381

TEL: +822 6902 8400 FAX: +822 6902 8401 www.osteonic.com

CS - E - 2009