




Optimus D PalaFix

Palatal Plate



Optimus D PalaFix



Palatal Plate



Features

- No risk of root damage, and no interruption of tooth movement
- Easy to perform with non-invasive treatment
- Teeth can be moved back without removal for protruded mouth in children during orthodontic procedures
- Fast and stable fixation compared to conventional methods
- Hook design prevents elastic chain droppage
- Having ample space for oral hygiene
- Dispenses with the need for incisions during plate placement
- Proper force is 400 ~ 500gf

Instrument

Item	Item Code
Screw - Ø2.0 x 8mm 	O19P08
Screw Driver Shaft - 35mm 	ODS-05

PAF - E - 1708